6th International conference on Advanced Nanomaterials,

Aveiro, Portugal, July 20-22 2015

Automatic equipment for hydrogen charge-discharge cycling and study of degradation and lifetime of hydrides

E. S. Ribeiro¹, F. J. A. Cardoso² and J. M. Gil¹,*

¹ CFisUC, Department of Physics, University of Coimbra, P-3004-516 Coimbra, Portugal ² GAII, Department of Physics, University of Coimbra, P-3004-516 Coimbra, Portugal *jmgil@fis.uc.pt

The study and characterization of new materials for use in hydrogen storage includes knowing their performance, stability, aging and/or degradation resulting from charge-discharge cycles. With this objective in mind, our group has developed an automated device whose purpose is to repeatedly undergo hydrogen charging and discharging in order to assess the longevity of the hydride. The equipment was tested with the LaNi₅-H system, where a significant degradation is observed over 1500 cycles.


Introduction

The characterization of the hydrogen absorption of new materials for hydrogen storage is typically performed with high accuracy Sieverts-type volumetric devices.


Studying the aging of hydrides requires taking a sample through many hydriding cycles [1], making conventional equipment inefficient for this task [2]. Here we present a dedicated device for the cycling over a significant number of cycles.

Pressure cycling system


The developed equipment is a volumetric system mounted with purchased elements from Swagelok Ltd, Keller Ltd. pressure transducer, and electrovalves from ASCO Numatics Ltd. The sample chamber is equipped with heating and cooling devices with temperature control in the range from 10°C to 230°C. The test material (LaNi₅, 99.9%) was purchased from Alfa Aesar.



A dedicated control and communication board developed for the automation this equipment. controlled experiment through a code developed in LabView with a graphical interface, including a PID thermal control for allowing the study different temperatures.







Acknowledgements


This work was partially supported by the Program COMPETE: PEst-C/FIS/UI0036/2011 strategic project and FCT through project PTDC/EME-MFE/103051/2008. E.S.R. wishes to acknowledge the Portuguese funding institution FCT - Fundação para a Ciência e Tecnologia for the Ph.D. Grant SFRH/BD/78733/2011.

PROGRAMA OPERACIONAL FACTORES DE COMPETITIVIDADI






Results




- a) Partial view of the pressure vs. time cycling of LaNi5;
- b) One hour hydriding/dehydriding cycle.

The shape of the absorption/desorption isotherm curves is also relevant to the degradation of hydriding properties of an intermetallic compound.



The hydrogen capacity loss as function of the number of pressure cycles for LaNi₅.



Conclusion

The testing equipment was successfully implemented, where we observe a loss of approximately 20% of the absolute capacity for LaNi₅-H system with 1500 cycles. For the study of new materials, we can use this equipment to simulate real situations of hydrogen storage tanks.

References

- [1] Lambert, S.W., et al. J. Alloys Compd, 187, 113 (1992).
- [2] G. Meyer, et. al. Rev. Sci. Instrum. 78, 023903 (2007).